

IN-OFFICE GEL - obj. kód 7343

Aplikace:

Nejprve naneste předbělicí lak, poté proveďte 3 aplikace BEAM In-Office gelu, kdy gel pokaždé nechte působit 15 minut.

Tento gel je určený pro použití s lampou BEAM.

Pro následnou domácí péči je určený BEAM WHITENING KIT.

Důležité informace k použití:

Přípravek není určen pro osoby mladší 18 let.

Výrobek je určen k prodeji pouze zubním lékařům / stomatologickým odborníkům.

Při každém cyklu použití musí být první aplikace provedena pouze zubním lékařem nebo pod jeho přímým dohledem, pokud je zajištěna odpovídající úroveň bezpečnosti. Následně je možné, aby spotřebitel výrobek převzal a dokončil finální cykly.

Není určeno pro těhotné nebo kojící ženy.

Zabraňte kontaktu s očima. Pokud se výrobek dostane do očí, okamžitě je vypláchněte vodou.

Uchovávejte na chladném a suchém místě při teplotě 5–25 °C, mimo přímé sluneční záření a mimo dosah dětí.

In-Office Whitening Gel (2,5 ml):

Složení: Aqua, Glycerin, Phthalimidoperoxycaproic Acid, Hydrogen Peroxide, Triethanolamine, Carbomer, Sodium Acryloyldimethyltaurate/VP Crosspolymer, Monosodium Citrate, PVM/MA Copolymer, Xanthan Gum, Sodium Saccharin, C12-15 Pareth-5.

Pre-Whitening Varnish (1 ml)

Složení: Glycerin, Aqua, Sodium Bicarbonate, Potassium Citrate, Hydroxyapatite, Glycolipids, Sodium Acryloyldimethyltaurate/VP Crosspolymer, Mentha Arvensis Leaf Oil, Stevioside, Mentha Piperita Oil, Potassium Hydroxide, Limonene.