

SMART | ONE FACEMASKS

Disposable 3-layer type IIR Face Masks with elastic band

CONFIGURATION OVERVIEW

- Designed for medical use
- Type IIR Face Mask
- Soft and comfortable
- Water repellant
- Fiberglass-free
- Latex-free
- PVC and Phthalates free
- Good performance of anti-micro-dust, anti-spray, anti-virus



Medical Device of Class I according to EU REG MDR 2017/745

CONFORM TO: EN ISO 10993-1:2021, EN ISO 10993-5:2009, EN ISO 10993-10:2010, ISO 15223-1:2021/Amd. 1:2025, EN ISO 14971:2020, EN ISO 20417:2021, EN 14683:2025

TECHNICAL SPECIFICATIONS:

- 3 layers = 2 polypropylene layers + 1 Melt blown layer (middle)
- Outer Layer: 35 g/m², PP, Coloured
- Middle Layer: 25 g/m², Purpose-made Meltblown, White
- Inner Layer: 25 g/m², PP nonwoven material, White
- 2 Holding Elastics: urethane elastic fiber
- Nose barret: Polyethylene with double iron wire bar

Dimensions of the mask (closed)	Dimensions of the mask (opened)	Length of nose barrette	Force at break for earloop
17.5 x 9.5 cm (± 0.5 cm)	17,5 x 15,5 cm (± 0.5 cm)	13 cm ±1 cm	1 Kg.

SHELF LIFE : 3 years

ARTICLE NUMBERS AND COLOUR VARIANTS:

FAMAA24-AKZ-BL-K02 = blue
 FAMAA24-AKZ-BP-K02 = baby pink
 FAMAA24-AKZ-FG-K02 = fresh green
 FAMAA24-AKZ-FU-K02 = fuchsia
 FAMAA24-AKZ-GR-K02 = green
 FAMAA24-AKZ-HO-K02 = orange
 FAMAA24-AKZ-LI-K02 = lilac
 FAMAA24-AKZ-MB-K02 = med blue
 FAMAA24-AKZ-PI-K02 = pink
 FAMAA24-AKZ-TL-K02 = laguna deep blue
 FAMAA24-AKZ-WH-K02 = white
 FAMAA24-AKZ-YE-K02 = yellow



Akzenta International SA
 Via Giuseppe Motta, 24
 CH-6830 Chiasso

tel. +41 (0) 91 921 14 92
 fax. +41 (0) 91 682 85 20

www.akzenta.com
 e-mail: business@akzenta.com

PACKAGING

Pieces/dispenser: 50
Dispenser/carton: 40
Carton size: 52 x 38 x 34 cm
Carton weight: 8,5 Kg
Cartons/layer : 4
Layers/pallet : 6
Cartons/pallet : 24
Pallet size : 80x120x219 cm
Pallet weight: 224 Kg

INSTRUCTIONS FOR USE:

1. Hold the mask and face the outer surface outward.
2. With the outer surface still facing away from you, place the ear loops comfortably around your ears.
3. Adjust the noseband to fit securely over the bridge of your nose, ensuring it conforms to the shape of your face.
4. While holding the noseband in place, pull the bottom edge of the mask under your chin, making sure it fully covers your nose, mouth, and chin.



Rev. 1 12-08-2025



Akzenta International SA
Via Giuseppe Motta, 24
CH-6830 Chiasso

tel. +41 (0) 91 921 14 92
fax. +41 (0) 91 682 85 20

www.akzenta.com
e-mail: business@akzenta.com